Clinical Barriers for Light Therapy Pads

Lightweight high clarity disposable bags or plastic wrap are placed over the limb or pad

Bag Source: ULINE.com

1 MIL weight is light, comfortable and economical: about 4- 6 cents a bag depending on the size. They are sold in boxes of 1000 and are made of 100% virgin, high clarity polyethylene film. Many more sizes and weights are available at www.uline.com.

2 MIL is also available if you prefer a heavier bag



Bags can be used in 2 ways:

- Place the limb in the bag and the pad over top as shown
- For all other uses: put the pad in a bag and place over the therapy location
- See page 2 for cap barriers for head therapy
- Plastic stretch wrap is also useful as a barrier, see page 2

 Boot/122 Pad For complete coverage over the leg and foot, place the bag on leg first, then strap the pad over top

10 x 24" 1 Mil Poly Bags S-6309

 Body/264 and Local/132 Pads, it can also cover a foot or hand. Cut shorter to use for 132 if desired, width is good

8 x 18" 1 Mil Poly Bags S-6307

Local/132 Pads, also use to cover hand

8 x 12" 1 Mil Poly Bags S-3641

Painstopper/180 Pad

6 x 18" 1 Mil Poly Bags S-10870

Plastic Cap Source: Amazon.com - Head Therapy Barrier

• Individually Wrapped clear shower cap for head therapy barrier

Zaveti Disposable Shower Caps - 200Pcs Waterproof Clear

Stretch Wrap Clear Plastic Barrier

This is the best barrier for the Facemask/104 and can be placed to allow any method of Velcro attachment

Clear plastic wrap may also be useful for wound healing applications

Stretch wrap can be used over therapy pads or therapy site

- PLASTIC WRAP can be placed **over the pad** much like wrapping a present.
- PLASTIC WRAP can be placed over the body part or area



A large roll from Costco is easy to handle. This works well for wound healing and skin applications. Use the optional cutter that comes in the box - the blade slides to cut the film.